



# actril<sup>®</sup> Cold Sterilant

## For Broad Spectrum Disinfection

- Sterilant
- Sporicidal
- Virucidal
- Tuberculocidal
- Fungicidal
- Hospital Level Disinfectant
- Cleaner-Sanitizer
- Ready-to-Use
- Non-Staining
- Fast-Acting

# Are You in Compliance?



FOR USE IN HEALTH CARE ENVIRONMENTS

Actril<sup>®</sup> Cold Sterilant is intended for use in Health Care Environments. When used as described in the enclosed spectrum disinfectant, hospital disinfectant, and germicidal spray disinfectant. As a sterilant Actril<sup>®</sup> Cold Sterilant is tuberculocidal, and will kill Hepatitis B and HIV-1 viruses. It is a violation of Federal Law to use this product in a health care environment.

**NET CONTENTS**  
**1 U.S. Gallons (3.78 Liters)**

Read Insert Carefully Before Using This Product

ACTIVE INGREDIENTS		
Hydrogen Peroxide	1.00%	EPA Reg. No. 52252-7
Peroxyacetic Acid	.08%	
Inert Ingredients	98.92%	EPA Est. No. 52252-MN-01
<b>Total</b>	<b>100.00%</b>	

**KEEP OUT OF REACH OF CHILDREN**  
**DANGER - PELIGRO**

**PRECAUTIONARY STATEMENTS**  
Hazard to Humans and Domestic Animals.

**DANGER - PELIGRO**

**Corrosive. Causes irreversible eye damage.** Harmful if absorbed through skin. Do not get in eyes, on skin, or on clothing. Avoid contact with skin. Wash hands before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash hands before reuse. Caution should be used when applying indoors because pets may be at risk.

Personal Protective Equipment  
Handlers must wear: Goggles  
Environmental Hazard  
This pesticide is toxic to fish and aquatic life.  
Physical and Chemical Hazards  
This product contains a flammable liquid.

Manufacturer  
**MINNTECH**  
A Corning Company  
Minntech Corporation  
14802 28th Avenue  
Minneapolis, MN 55425  
Phone: 763-833-1000  
Toll Free: 1-800-451-1000  
Fax: 763-833-1001



U.S. Pat. No. 5,656,302 and foreign patents

Distributed By:



**MINNTECH<sup>®</sup>**  
**renal systems**

A Minntech Corporation Business Group

# Actril<sup>®</sup> Cold Sterilant

# actril<sup>®</sup>

## Overview

Actril<sup>®</sup> Cold Sterilant is a stabilized mixture of peracetic acid, hydrogen peroxide and acetic acid, which leaves no toxic residue. For two decades, Actril<sup>®</sup> Cold Sterilant continues to give clinics everything they could ask for in a surface disinfectant.

Actril is a powerful broad spectrum disinfectant that is safe to use, environmentally friendly, leaves no toxic residue and is thoroughly proven in hundreds of clinics.



**Actril<sup>®</sup> Cold Sterilant System** (4 U.S. Gallons [3.78 liters] per case). Includes 100 each Residual and Indicator Strips  
**P/N 78270-000**

**Actril<sup>®</sup> Cold Sterilant** (4X1 U.S. Gallons [3.78liters] per case)  
**P/N 78337-000**

**Actril<sup>®</sup> Residual Test Strips** (100 per vial)  
**P/N 78258-000**

**Actril<sup>®</sup> Indicator Test Strips** (500 per pack)  
**P/N 78259-000**

## Accuracy and Validation

• **Surface Disinfection** – Actril<sup>®</sup> Cold Sterilant is a ready-to-use surface disinfectant, which is extremely effective, even in high bio-load environments. Unlike bleach, which quickly “deactivates” in the presence of organic material such as blood, Actril<sup>®</sup> Cold Sterilant chemistry continues to form peracetic acid even as it works to disinfect. Actril Cold Sterilant does not leave a toxic residue

• **Machine Disinfection** – Actril<sup>®</sup> Cold Sterilant has been successfully used for two decades for the disinfection of dialysis machines. Today, Actril<sup>®</sup> Cold Sterilant continues to be the ideal chemical for machine disinfection: effective against even the most resistant organisms, biodegradable, ready-to-use and safe. Several dialysis machine manufacturers have recognized the benefits of using Actril<sup>®</sup> Cold Sterilant and provide instructions for its use in their equipment. Additional information is available by contacting Minntech Renal Systems.

• **Reduces Scale Deposits** – Actril<sup>®</sup> Cold Sterilant contains acetic acid. Acetic acid has been used for decades to reduce or eliminate “scale” build-up (precipitated calcium and magnesium carbonate) when used on a regular basis. Actril Cold Sterilant also breaks down after use, leaving little residual that could interfere with safe machine operation.

## About Minntech

Minntech Renal Systems has been providing WORKING SOLUTIONS to the hemodialysis community since 1972. The Minntech Renal Systems lines of reprocessing systems, sterilants, disinfectants and dialysate concentrates have been used in thousands of clinics and tens of millions of treatments.

Today, Minntech Renal Systems continues to provide dialysis clinics with the highest quality products and services to ensure that clinics have the WORKING SOLUTIONS necessary for optimal patient care.

**ALSO AVAILABLE**  
**Actril<sup>®</sup> Cold Sterilant Machine Tag** (1 each)  
**P/N 67197-003**



Distributed By:



Minntech Corporation  
14605 28th Avenue North  
Minneapolis MN 55447  
Tel: +1-763-553-3300  
Fax: +1-763-553-3387  
www.minntech.com

Minntech BV  
Sourethweg 11  
6422 PC Heerlen  
The Netherlands  
Tel: +31-45-5-471-471  
Fax: +31-45-5-429-695

Minntech Asia/Pacific Pte Ltd  
438 Alexandra Road  
#03-01 Alexandra Point  
Singapore 119958  
Tel: +65-6227-9698  
Fax: +65-6225-6848

Minntech China, Rep. Office  
Room 1508B, Avic Plaza  
No. 10 Yi Dongsanhuan Zhonglu  
Chaoyang District,  
Beijing 100022 China  
Tel: +86-10-6567-8446  
Fax: +86-10-6567-8445

For More Information Contact Minntech at  
[www.minntech.com/renal](http://www.minntech.com/renal) or Call 800-328-3340