

SANIDOSE[®] TEST STRIPS

INSTRUCTIONS FOR STORAGE & USE

INTRODUCTION

The SANIDOSE Test Strips are easy to use and are conveniently packaged 100 test strips per container. Each strip has a "reaction zone" pad located at one end of the strip. When stored and used properly the test strip pad can indicate acceptable SANIDOSE levels.

STORAGE

- ✓ STORE CONTAINERS IN COOL AND DRY PLACE WITHOUT REFRIGERATING.
- ✓ FAILURE TO STORE TEST STRIPS PROPERLY WILL RESULT IN INACCURATE READINGS.

INSTRUCTIONS FOR USE

1. Collect a sample of the rinse solution.
2. Remove a test strip from the tube and immediately replace the lid.
3. Dip the test strip into the solution to be tested for one second, such that the reaction zone is properly wetted.
4. Remove the test strip, gently shake off excess liquid, and compare the reaction zone with the color scale between 5 and 10 seconds after exposure.

NOTE:

- If test strip indicates dark blue to brown or green to brown, the concentrations are acceptable.
- Read SANIDOSE Package Insert instructions for additional information.

Endoscope reprocessing the way it should be™

MEDIVATORS[®]

Manufactured in the USA by:



MINNTECH CORPORATION
Minneapolis, MN USA
Toll Free: (800) 444-4729

MINNTECH BV 
The Netherlands
Tel: +31 045 5 471471