

MEDIVATORS®  
**Rapicide® OPA/28**  
HIGH-LEVEL DISINFECTANT

## RAPICIDE OPA/28 Test Strips Quality Control Procedure


**\*Refer to the Directions for Use for more information**

### 1. Preparation of Control Solutions

To prepare positive and negative control solutions for testing, obtain a container of RAPICIDE OPA/28 High-Level Disinfect solution. Verify that the solution has not met or exceeded the labeled expiration date. For the positive control, the full strength solution should be used. To prepare a negative control, dilute one part of full strength RAPICIDE OPA/28 High-Level Disinfectant solution with one part of water at 95°F (35°C). Label each control solution appropriately.

### 2. Testing Procedure

Following the directions for use, dip one indicator strips in each of the above freshly prepared solutions. The strip dipped in the full strength positive control solution should exhibit a complete green color on the indicator pad when read at 90 seconds. The strip dipped in the diluted negative control should display a blue color when read at 90 seconds. Refer to the visual standard on the test strip bottle for interpretation of the results. This visual is also displayed below.



### 3. Testing Frequency

It is recommended that the testing of positive and negative controls be performed on each newly opened test strip bottle of RAPICIDE OPA/28 Indicator test strips. After this initial testing, it is recommended that testing of freshly prepared positive and negative controls be performed on a regular basis as established by your own quality control procedures and program. This testing program will serve to minimize errors between different users, use of outdated materials or product that has been improperly stored or handled.

### 4. Unsatisfactory QC Test performance

If the results obtained from using the positive and negative controls indicate the test strip is not functioning properly, discard the remaining strips. **DO NOT USE.** For customer support, contact MEDIVATORS at 1-800-444-4729.

MEDIVATORS® and RAPICIDE® is a registered trademark of Medivators Inc.

To order, call MEDIVATORS directly at: 1-800-444-4729, option 1 • [www.medivators.com](http://www.medivators.com)

