


Using MEDIVATORS® ADVANTAGE® Automated Endoscope Reprocessors

Locate the endoscope port and follow the instructions for attaching the corresponding Endoscope Connector. All connectors must be attached.

ADVANTAGE HOOKUP DIAGRAM



Each ADVANTAGE Hookup is comprised of multiple ports and connector assemblies. Each port is labeled with a Port ID number (circled). This same number also identifies the corresponding Connector Assembly.

CONNECTOR ASSEMBLY IDENTIFICATION CHART

Port ID	Port Name	Connector Assembly Description	Order Number
①	LEAK TEST	Leak Test Connector Assembly	H28511-1
②	JET/AUX	Air/Water Connector Assembly	H28511-26
③	BIOPSY	Biopsy Connector Assembly	H28511-3
④	AIR	Air Probe Connector Assembly	H28511-4
⑥	WATER	Air/Water Connector Assembly	H28511-26
⑦	SUCTION	Suction Connector Assembly	H28511-7
Ⓢ	NA	Air/Water & Suction Port Separator	78400-930

If a Connector Assembly needs replacement, locate the port number/name and use the corresponding order number found in the chart above. Connector Assemblies are specific to individual hookups.

If you have any questions regarding this hookup, in the U.S., please contact MEDIVATORS Technical Support at 1-800-444-4729. Outside the U.S., please contact your MEDIVATORS Distributor.

Air/Water and Suction Port Separator

Align the separator with piston ports. Ensure that the release button is pushed in. Push connector onto piston ports until the release button moves out and clicks.

Biopsy Connector

Align the biopsy connector disc over the biopsy inlet and push down to compress spring and slide the frame forward to fully engage.

Suction Connector

Align silicone tubing with suction barb and push tubing onto barb until fully engaged.

Leak Test Connector

Align the slot on connector with the pin on the scope. Push on and turn clockwise to engage. To remove leak test connector, turn connector counter-clockwise until slot is aligned with pin and pull off. Do not pull until slot is aligned with pin.

Air Probe Connector
Align with Air Probe, push and slide until fully engaged.

Air/Water Connector

Align the connector openings with the air/water ports of the endoscope as shown. Ensure side latch is in. Push connector onto ports until latch moves out and clicks. To remove, push latch to release.

DESIGNATED HOOKUP FOR THE FOLLOWING OLYMPUS® ENDOSCOPES*:

CF-140	CF-240I	CF-H260AZI	CF-Q240AL	GIF-240	GIF-PQ260L	GIF-Q165	GIF-SP240	GIF-XQ140	PCF-140L	PCF-240S	SIF-Q180	CF-Q160ZL	GIF-NP190N
CF-140I	CF-240L	CF-H260AZL	CF-Q240I	GIF-FQ260Z	GIF-Q140	GIF-Q180	GIF-NP150N	GIF-XQ160	PCF-160AI	PCF-P240AI	SIF-Q240	CF-Q260AI	GIF-XP170N
CF-140L	CF-40I	CF-Q140I	CF-Q240ZI	GIF-H180	GIF-Q145	GIF-Q240	GIF-XP160	GIF-XQ240	PCF-160AL	PCF-P240AL	SIF-Q245	CF-Q260AL	GIF-H170
CF-140S	CF-40L	CF-Q140L	CF-Q240ZL	GIF-H260	GIF-Q150	GIF-Q240X	GIF-XP180N	GIF-XQ260	PCF-240AI	PCF-Q260AI	SIF-Q260	CF-Q260DI	PCF-PH190I
CF-240AL	CF-H260AI	CF-Q150L	GIF-140	GIF-N180	GIF-Q160	GIF-Q240Z	GIF-XP240	GIF-XQ40	PCF-240AL	PCF-Q260AL	TJF-Q180V	CF-Q260DL	PCF-PH190L
CF-240AL	CF-H260AL	CF-Q160DL	GIF-160	GIF-P140	GIF-Q160Y9	GIF-Q260	GIF-XP260	PCF-140	PCF-240I	PCF-S	CF-Q160Z	GIF-H185N	
CF-240DL	CF-H260AXL	CF-Q240AI	GIF-XP290N	GIF-Q160Y9	GIF-Q160Z	GIF-Q40	GIF-XP260NS	PCF-140I	PCF-240L	SIF-Q140	CF-Q160ZI	GIF-H290	

* This list of corresponding endoscope models is subject to change without notice. For a full list, please go to www.medivators.com/hookuplookup




Ensure proper alignment of all connectors with the endoscope mating parts prior to attachment.

Representative endoscope for illustrative purposes only. For exact port locations, refer to endoscope manufacturer's instruction guide.

HOOKUP CONNECTOR BLOCK LABEL

The Hookup Connector Block Label contains useful information for quick reference.



Handle Locking Instructions - insert block onto basin channel manifold and move handle to the left to lock. Move to the right to unlock.

Lot Number - identifies date of manufacture.

Hookup Part Number - identifies hookup.

Hookup Barcode - barcode of part number.

Port Name - identifier of port location in ADVANTAGE system. Error logs identify issues using these names. Also identifier of Connector Assembly replacements for older model hookups.

Port ID Number - identifier for Connector Assembly replacements. Not found on older model hookups.

Alternative Channel Connector Names - connector names vary depending on specific endoscope channel configuration.

WARNINGS & CAUTIONS



This hookup is intended for connection only with the specific endoscope models identified in the appropriate MEDIVATORS online hookup guide. Use of this hookup with endoscopes other than those specified by the hookup application guide may result in inadequate HLD.

Users operating MEDIVATORS AERs and hookups must be trained and competent in the understanding of endoscope channel systems.

Prior to reprocessing, users must verify that the hookup used contains connections for all appropriate channels that require a separate reprocessing connection according to the hookup instructions

All connectors and adapters must remain firmly attached and unrestricted for the entire disinfection cycle to ensure adequate HLD. Users must inspect the endoscope and hookup combination to verify proper connection and flow.

Modifications or repairs to hookups which do not correspond to the manufacturer's specifications may result in inadequate HLD and/or damage to endoscopes.

Failure to properly clean and prepare endoscope for immersion prior to hookup connection may result in inadequate HLD and/or damage to the endoscope.

Endoscopes must be inspected for damage and verified to be in proper working order prior to hookup connection. Connection of hookups to damaged endoscopes or endoscopes with obstructed/restricted channels could result in inadequate HLD.

Manual leak testing must be carried out in accordance with manufacturer's recommendations prior to reprocessing.

Install all necessary waterproof caps, plugs and cleaning adapters prior to immersion.

Do NOT autoclave MEDIVATORS hookups.

MEDIVATORS makes no claim on high-level disinfection efficacy when these reprocessing instructions are not followed, or when this hookup kit is applied to endoscopes other than those specified in the online interactive hookup guide.

Avoid biological contamination. Always wear personal protective equipment when handling endoscopes.



Manufactured in the USA by:

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Toll Free: +1.800.444.4729



Medivators BV
Sourethweg 11
6422PC Heerlen
The Netherlands
Tel: +31.45.5.471.471



Cantel Medical Asia/Pacific Pte. Ltd.
1A International Business Park
#05-01 Singapore 609933
Tel: +65.6227.9698

Cantel Medical Devices (China) Co. Ltd.
Unit 804-805, Innov Tower Block A,
Hongmei Road, Xuhui 200233 Shanghai
Tel: +86 21 60161380
Fax: +86 21 61210913

www.medivators.com