

	Actions	Facility Installs	Medivators
Site Inspection and Evaluation	Assess location for space, electrical, water, drain, electrical and ventilation	Assist on request	Assist on request
	Design/reconfigure area for work flow optimization	Assist on request	Assist on request
	Evaluate water for hardness, particulates	Assist on request	Assist on request
Site Preparation	Provide electrical supply	●	●
	Provide water supply	●	●
	Connect water pre-filtration	●	●
	Provide drain connection	●	●
	Provide adequate ventilation	●	●
Equipment Setup	Uncrate Equipment	●	●
	Place in designated space	●	●
	Level	●	●
	Connect to power/water/air/drain	●	●
	Attach VMS (if supplied)	●	●
	Install data management system (if supplied)	●	●
	Load disinfectant	●	●
Performance Certification	Verify proper installation	Included	●
	Re-program for alternate LCG (if applicable)	Included	●
	Perform test cycles	Included	●
	Activate warranty	Included	●
	Conduct "quick-start" operator training	Included	●
	Review "quick start" reference materials	Included	●
Operator Training	Follow-up staff training on AER operation/troubleshooting	Included	●
	Additional comprehensive staff training	●	●

■ Facility Responsibility
 ■ Medivators provides at no charge
 ■ Medivators provides for a fee

For Technical Support, call 1-800-444-4729



Customer support,
service & training

Endoscope reprocessing the way it should be

