

CENTRISOL® “100 Dextrose” Concentrate Formulations (Part A)

These low dextrose (100 mg/dL) formulas were developed upon request by our customers. They are made according to the same quality standards as all of our formulations. Contact your Minntech Renal Systems’ representative regarding availability in your area.



Part A: SB-300 Series Liquid Acid Proportioning Components

Formula	Sodium mEq/L	Calcium mEq/L	Potassium mEq/L	Magnesium mEq/L	Chloride mEq/L	Acetate mEq/L	Dextrose mg/dL	Bicarbonate mEq/L	Packaging Cases	Drums
SB-330	137	0.00	0.00	0.75	100.75	4.00	100	33.00	X	
SB-351	137	2.00	2.00	0.75	104.75	4.00	100	33.00	X	X
SB-352	137	2.50	2.00	0.75	105.25	4.00	100	33.00	X	X
Low Calcium High Acetic Formulas										
SB-342	137	2.25	1.00	0.75	104.00	5.00	100	32.00	X	X
SB-343	137	2.25	2.00	0.75	105.00	5.00	100	32.00		X
SB-344	137	2.25	3.00	0.75	106.00	5.00	100	32.00	X	X

Quality, Everyday.



**Concentrate Formulation Questions?
Call Clinical Services 1-800-328-3345**



MINNTECH®
renal systems

Terms: F.O.B. Minneapolis, MN. Effective October 1, 2010. Call your Minntech Representative for current pricing information.