

CVVHD Clearances – HF 700

Clearances versus inlet dialysate flow rate

(Continuous venovenous hemodialysis) Saline, T=37°C

	QB 100 (ml/min)				QB 200 (ml/min)			
QD (ml/min)	17	42	67	133	17	42	67	133
Urea (mg/dl) ($\pm 10\%$)	19	35	51	70	20	38	63	86
Vitamin B-12 (mg/dl) ($\pm 10\%$)	16	28	36	45	18	30	40	52
Cytochrome-C (mg/dl) ($\pm 10\%$)	15	21	24	26	15	20	24	25