



Dimensions

9.5 H x 6 W x 8 D (inches)
 (24.1 H x 15.2 W x 20.3 D (cm))

Weight

10 lbs (4.5 kg)

Electrical Requirements

110 – 230 V, 50/60 Hz

Electrical Safety Certifications

UL Standard 61010-1
 CANADA/USA C22.2 No. 61010-1
 ETL and cETL listed

Fluid Delivery Rate

1 liter/minute

Flush Time

0 – 99 minutes

Timer Control

Digital display with touchpad controls

System Pressure

Built-in safety bypass

Cycle Completion Indicators

Audible and visual alerts

Unit Placement

Counter-top or wall-mounting (bracket included)

Warranty

One year

Scope Buddy Ordering Information

Part number	Description
EFA-US-G	North America
EFA-IN-G	International

Scope Buddy Accessories

Tubing Kit

Complete set of connection tubings
 Part Number 78399-669

Intercept® Detergent

Detergent solution used during endoscope flushing.
 Part Number ML02-0106

Graduated Cylinder

250ml cylinder used during quality assurance test
 Part Number 75227-679



Minntech Corporation
 14605 28th Avenue North
 Minneapolis MN 55447
 Tel: +1-763-553-3300
 Fax: +1-763-553-3387
 www.minntech.com

Minntech Asia/Pacific Pte Ltd
 No.1A International Business Park
 #05-01 Singapore 609933
 Tel: +65-6227-9698
 Fax: +65-6225-6848

Minntech BV
 Sourethweg 11
 6422 PC Heerlen
 The Netherlands
 Tel: +31-45-5-471-471
 Fax: +31-45-5-429-695

Minntech China, Rep. Office
 Room 1508B, Avic Plaza
 No. 10 Yi Dongsanhuan Zhonglu
 Chaoyang District,
 Beijing 100022 China
 Tel: +86-10-6567-8446
 Fax: +86-10-6567-8445

For additional information please contact your Medivators representative
 or visit our website at www.minntech.com/medivators

Eliminate manual flushing of endoscopes



Endoscope reprocessing the way it should be

MEDIVATORS

SCOPE BUDDY®

Consistent, hands-free flushing

Reduces the risk of repetitive motion injuries

- Replaces repetitive and labor-intensive syringing of channels
- Hands-free operation reduces muscle fatigue
- Push button operation
- Provides constant fluid delivery

Enhances compliance with endoscope cleaning protocols

- Continuous, timed fluid delivery ensures consistent, repeatable endoscope flushing
- Audible and visual alerts upon completion
- Compatible with all flexible endoscopes including ultrasound and those with elevator wire channels
- Works with your existing sink; compact size accommodates countertop placement of wall mounting
- Improves manual cleaning throughput

Shortens manual cleaning time

- Thoroughly flushes most endoscopes in 35 seconds
- Flushes multiple endoscope channels simultaneously
- Quickly and quietly delivers detergent and rinse solutions through endoscope channels

Combine Intercept® detergent with Scope Buddy for...

- Effective endoscope bacterial biofilm penetration and removal
- Superior removal of biological and organic soils
- Assured endoscope compatibility



ng of endoscopes

Visual display of endoscope flushing time

User-selectable time setting

One-touch start/stop operation

Simple and quick tubing connections

