

Facility Name							
Address							
City		State		Zip Code		Country	
Installation Date				Installed by			
Sales Representative							

This checklist will be utilized to ensure the customer installation requirements are met for the DSD Data Management Software .

INSTALLATION CHECKLIST STANDARD INSTALLATON	
<input type="checkbox"/>	Customer supplied Dedicated PC (or PC's depending on number of machines connected)
<input type="checkbox"/>	PC is running Windows XP
<input type="checkbox"/>	PC has Windows 3.1 Installer (or newer)
<input type="checkbox"/>	Microsoft.NET Framework version 2.0.50727 (or newer)
<input type="checkbox"/>	PC station needs to be less then 40 – 45 ft from the DSD (should be a dedicated PC)
<input type="checkbox"/>	PC has 1 serial port DB9 per DSD being connected (and any additional hardware for more machines)
<input type="checkbox"/>	NIC Card (If Networked)
<input type="checkbox"/>	A CD-ROM Drive (to load software)
<input type="checkbox"/>	If PC is located beyond 40 – 45 ft. Completed connectivity from PC to DSD
NOTES	

Facility Representative:			
Title:			
Signature:		Date:	

Provide a copy to the Facility Representative and Fax the completed form to:
Medivators Reprocessing Systems | Technical Support | FAX (866) 553-9444
 Minntech Corporation 14605 28th Avenue North, Minneapolis, MN 55447 USA