



PRODUCT BULLETIN

<i>Distribution:</i>	
International Distributors	✓
International Management	✓
US Clinical Specialists	✓
US Customer Service	✓
US Field Sales	✓
US Technical and Field Service	✓

Date: January 18, 2011

RE: Preventative Maintenance Programs for Medivators DSD, SSD and Advantage Series Reprocessors

As with any sophisticated electro-mechanical device that must operate and perform with stringent performance specifications, endoscope reprocessors (AERs) should receive regular maintenance checks to ensure proper and safe operation.

The components used within these devices have been carefully selected to meet the operating requirements of the environment in which they perform. Many of the components do however, have a finite life beyond which they will deteriorate, especially when used with elevated temperatures.

Many parts perform metering functions, can wear out and will cause incorrect amounts of fluids to be dosed in various phases of a disinfection cycle. Worn parts can also increase the amount of disinfectant dilution, which will shorten the re-use life of the high level disinfectant in use.

AERs should also be checked for cycle functionality, software revision, fume containment and many other parameters that ensure proper operation. At a minimum, a Preventative Maintenance service is recommended annually. High volume endoscopy facilities may require a more frequent check.

Please contact your local Medivators Sales consultant, or the Medivators Customer Support Group at 1-800-444-4729 for full details of available maintenance agreements.

Endoscope reprocessing the way it should be

MEDIVATORS®

Minntech Corporation, 14605 28th Avenue North, Minneapolis, MN 55447-4822

763.553.3300, fax 763.553.3387, www.minntech.com