



# Technical Bulletin

<i>Distribution:</i>	
Olympus America Inc.	<input type="checkbox"/>
International distributors	<input type="checkbox"/>
MediVators – Field Clinical Specialists	<input type="checkbox"/>
MediVators – Technical Service	<input type="checkbox"/>

**Date:** October 16, 2001

**RE:** Pentax Ultrasound Endoscopes

## Distributor/Customer

Certain Pentax Ultrasound endoscope models have a scanning unit connector that cannot be immersed in any form of solution. This connector is **not** waterproof.

To ensure that these endoscopes can be properly reprocessed in the MediVators DSD/SSD range of reprocessors, MediVators has produced a support cradle to enable this connector to be positioned out of contact of solutions being used during a disinfection cycle.

The part number for this kit is MK01-0053 and consists of the support bridge and slotted floating lid.

If you have any questions relating to the contents of this Technical Bulletin, please contact your representative or MediVators Inc.

Nick Bromiley  
Senior Technical Specialist