

<i>Distribution – Internal Use Only</i>	
International Distributors	✓
International Management	✓
US Clinical Specialists	✓
US Customer Service	✓
US Field Sales	✓
US Technical and Field Service	✓

Date: November 1st, 2009

RE: CER/DSD hook-up changes

Over the past year, numerous connectors used on CER and DSD hook-ups have been redesigned for one or more of the following reasons:



- Ease of use improvement
- Increased use life
- Vendor change






Implementation of these connector changes into corresponding hook-ups has been prioritized and completed for the highest volume CER/DSD hook-ups. The remaining affected CER/DSD hook-ups will be updated with the most recent connector versions by the end of December 2009. All existing inventory on-hand in our warehouse will also be reworked. A list of connectors and the reason for their redesign follows, along with pictures of the connectors.

During the process, Directions For Use will be updated for all CER/DSD hook-ups.



Upon release of an updated hook-up, a Product Bulletin will be distributed, citing the connector changes, and a copy of the new DFU will be provided.

If you have any questions related to this bulletin, please contact Medivators Technical Support Department at 1-800-444-4729. Outside of the United States, please contact your local Medivators distributor.

Item (<i>newest connector on right</i>)	Name/Description	Reason for Redesign
	Olympus Double Biopsy Channel Inlet Cap	Vendor Change
	Single Biopsy Channel Inlet Cap	Vendor Change
	Olympus Molded 10/20 Series Gas Cap	Vendor Change
	Bronchoscope Biopsy Inlet Connector	Vendor Change
	Air Probe, Olympus	Increase use life

Item (<i>newest connector on right</i>)	Name/Description	Reason for Redesign
	Suction Connector	Develop a universal connector for Olympus, Fujinon and Pentax; improve ease of connection to scope
	Auxiliary Water Channel Connector for Olympus 130 series and below	Vendor Change
	Universal Luer	Increase use life
	Biopsy Channel Adapter Connector	Vendor Change
	Elevator Wire Channel all Olympus and Auxiliary Water Channel Connector, Olympus 140 series and above	Improve ease of connection to scope

Item (<i>newest connector on right</i>)	Name/Description	Reason for Redesign
	Suction Adapter, Olympus	Vendor Change
	Leak Test Adapter Fujinon	Re-designed to fit all series scopes, including 530 series
	Air Probe, Fujinon	Improve security of attachment to Fujinon air probe connection
	Air/Water/Gas (AWG) Plug for Olympus 130 series and above	Vendor change
	Air/Water Connector Plug for Olympus 140 series and above	Increase use life

Item (<i>newest connector on right</i>)	Name/Description	Reason for Redesign
	Air/Water Channel Plug for Olympus	Vendor change
	Air/Water Channel Plug, Olympus Ultrasound	Vendor change