

# ADDENDUM

## **Subject: Commercial Solutions® Formula 409® Heavy-Duty Degreaser**

### **Notice:**

The Clorox® Company discontinued the Formula 409 NQF For Institutional Use. Therefore, Commercial Solutions Formula 409 Heavy-Duty Degreaser should be used to perform the Formula 409 Service Cleaning Procedure.

### **WARNING:**

Only the Commercial Solutions Formula 409 Heavy-Duty Degreaser should be used to perform the Formula 409 Service Cleaning Procedure. Other cleaners or degreasers may damage the Renatron station. If your facility still has Formula 409 NQF For Institutional Use in stock, it is still an acceptable chemical to use in the Renatron® station.

If you have any questions, please call MEDIVATORS Technical Service at 1-800-328-3324. If you would like to order the Commercial Solutions Formula 409 Heavy-Duty Degreaser, please call MEDIVATORS Customer Service at 1-800-328-3340.