

MEDIVATORS™ ADVANTAGE

Hookup Rack

The MEDIVATORS ADVANTAGE Hookup Rack stores hookup blocks used in the ADVANTAGE® Endoscope Reprocessor. The graduated design keeps hookup blocks organized for quick retrieval and easier tube management. The hanging style improves the ability of the hookup tubing and connectors to dry and removes countertop clutter caused by hookups.

- Hanging for improved tube and connector drying
- Graduated stacking holds more hookups in an area and reduces tube tangling
- All hookup labels face forward for quick identification
- Storage row designed specifically for disinfection and soak-only blocks
- Convenient shelf holds channel separators and accessories
- Stainless steel

ADVANTAGE Hookup Rack Ordering Information

Model number	Description
78400-915	Holds 16 hookup blocks plus 4 disinfection/soak only blocks Dimensions: 22"W x 15.75" D x 24.5" H
78400-931	Holds 8 hookup blocks plus 2 disinfection/soak only blocks Dimensions: 11.25" W x 15.75" D x 24.5" H



MEDIVATORS™ HOOKUP RACK

The MEDIVATORS Hookup Rack is designed for DSD & CER automated reprocessor hookups, flushing aid tubing for use with SCOPE BUDDY® Flushing Aids and tubing sets compatible with VERISCAN® Leak Testers. This convenient storage rack organizes hookups and tubing sets and improves drying. Optional shelf style holds channel separators and other accessories.

- Fits DSD and CER reprocessor hookup tubing and tubing sets for SCOPE BUDDY Flushing Aids and VERISCAN Leak Testers
- 9 slots each hold up to 3 DSD hookups, 4 CER hookups, 3 tubes for SCOPE BUDDY Flushing Aid and 2 tubes for VERISCAN Leak Tester
- Slots designed for easy in and out

Hookup Rack Ordering Information

Model number	Description
72233-774	MEDIVATORS Hookup Rack Dimensions: 13"W x 2.5" D x 3" H



78400-894	MEDIVATORS Hookup Rack with Shelf Dimensions: 13.375" W x 6" D x 4" H
-----------	--



MEDIVATORS Inc. Headquarters
14605 28th Avenue North
Minneapolis, MN
55447-4822 USA
Tel: +1.800.444.4729
Fax: +1.763.553.3387

MEDIVATORS BV
Sourethweg 11
6422 PC Heerlen
The Netherlands
Tel: +31.45.5.471.471
Fax: +31.45.5.429.695

MEDIVATORS Asia/Pacific Pte. Ltd.
1 A International Business Park
#05-01 Singapore 609933
Tel: +65.6227.9698
Fax: +65.6225.6848

MEDIVATORS Beijing Representative Office
Room 708, 7th Floor Kaiheng Center, Block B
No. 2 Chaoyangmennei Street
Dongcheng District, Beijing China 100010
Tel: +8610.6567.8446
Fax: +8610.6567.8445

www.medivators.com

© 2012 MEDIVATORS Inc. P/N 50097-976/C