

MEDIVATORS®  
THE **ADVANTAGE**  
POLYP TRAP® COLLECTION CONTAINER



# Easy, In-Line Polyp Recovery

The MEDIVATORS® THE ADVANTAGE POLYP TRAP® and ECONOTRAP® Polyp Traps are designed for easy in-line installation to facilitate uninterrupted suction and eliminate the need to retrieve polyps from collected endoscopic effluent. The removable numbered chambers allow for accurate specimen identification and minimize the exposure to collected biomaterials.

## MEDIVATORS® THE ADVANTAGE POLYP TRAP® COLLECTION CONTAINER



- Reduce risk of lost colonic tissue
- Eliminate manual manipulation after retrieval
- Large numbered chambers for quick, easy transfer and accurate specimen identification
- Single-use to protect your patients from potential cross contamination

## MEDIVATORS® ECONOTRAP® POLYP TRAP



## MEDIVATORS® ECONOTRAP®<sup>2</sup> POLYP TRAP



### Ordering Information

Model Number	Description	QTY: 50/Box
100101	ECONOTRAP Single Chamber Polyp Trap	QTY: 50/Box
100102	ECONOTRAP2 Double Chamber Polyp Trap	QTY: 50/Box
100104	THE ADVANTAGE POLYP TRAP Four Chambered Polyp Trap	QTY: 50/Box

MEDIVATORS®, THE ADVANTAGE POLYP TRAP® and ECONOTRAP® are registered trademarks of Medivators Inc.



Medivators Inc.  
14605 28th Avenue North  
Minneapolis, MN  
55447-4822 USA  
Tel: +1.763.553.3300  
Fax: +1.763.553.3387

Medivators BV  
Sourethweg 11  
6422 PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471  
Fax: +31.45.5.429.595

Cantel Medical (UK) Limited  
Wulseley House  
Staffordshire Technology Park  
Beaconside, Stafford ST18 0GA  
United Kingdom  
Tel: +44 (0) 1785.782420  
Fax: +44 (0) 1785.782421

Cantel Medical Asia/Pacific Pte. Ltd.  
1A International Business Park  
#05-01 Singapore 609933  
Tel: +65.6227.9698  
Fax: +65.6225.6848

Medivators Inc. Beijing Representative Office  
Room 1801, Floor 18th, Tower A,  
Beijing Marriott Hotel, Office Building  
No. 7, Jianguomen South Avenue,  
Dongcheng District, Beijing 100005 China  
Tel: +86 10 65204039

[www.medivators.com](http://www.medivators.com)