

March 11, 2011

Re: Medivators AERs High Level Disinfection of the KARL STORZ Flexible GI Videoscopes

Dear Valued Customer:

KARL STORZ Endoscopy-America, Inc. shares the responsibility with our customers of providing the best reprocessing practices. KARL STORZ, in collaboration with Minntech Corp., has completed the high level disinfection validations to represent reprocessing of the KARL STORZ Flexible GI Videoscopes. The high level disinfection validations met the required 6-log reduction of the biological indicators, *Mycobacterium terrae* ATCC # 15755, using the representative Medivators AER systems.

The Medivators AER table matrix (attached) lists the qualified KARL STORZ GI Videoscopes, Medivators Automated Endoscope Reprocessors (AER) and specific Medivators hook-up connectors.

Note: The Dual Sided Disinfector (DSD), Single Sided Disinfector (SSD), and DSD Edge Hook-up Connectors require the KARL STORZ Washer/Disinfector Connector Set 13991CK.

Should you have any further questions regarding this subject or require additional assistance, please don't hesitate to KARL STORZ Technical Support at (800) 421-0837. You may also contact Minntech Corp. Technical Support at (800) 444-4729. Minntech Corp. is the manufacturer of Medivators AERs.

Sincerely,



Tony Rosales
Sterilization Manager
Regulatory Affairs Department

Ref.: 2011-012-AR

**Medivators Automated Endoscope Reprocessors
KARL STORZ Flexible GI Videoscopes**

KARL STORZ Flexible GI Videoscopes		Medivators AERs and Hook-up Connectors					
Description	Model No.	CER-1 (MV-1)	CER-2 (MV-2)	DSD / SSD	DSD Edge	Advantage	Advantage Plus
<i>Single Channel</i>		CLM-110-HU0080 (x4) and KARL STORZ Connector Set 13991CK	CLM-110-HU0080 (x3), CLM-110-HU0136, and KARL STORZ Connector Set 13991CK	DSD-110 -HU0163 and KARL STORZ Connector Set 13991CK		2-8-540HAN	
Gastrosopes	13800 -NKS/-PKS						
	13801 -NKS/-PKS						
	13802 -NKS/-PKS						
	13803 -NKS/-PKS						
	13804 -NKS/-PKS						
	13805 -NKS/-PKS						
	13808 -NKS/-PKS						
	13820 -NKS/-PKS						
Duodenoscope	13821 -NKS/-PKS						
	13880 -NKS/-PKS						
	13881 -NKS/-PKS						
	13882 -NKS/-PKS						
Sigmoidoscopes	13885 -NKS/-PKS						
	13900 -NKS/-PKS						
Colonoscopes	13901 -NKS/-PKS						
	13902 -NKS/-PKS						
	13903 -NKS/-PKS						
	13904 -NKS/-PKS						
	13905 -NKS/-PKS						
	13906 -NKS/-PKS						
	13907 -NKS/-PKS						
13924 -NKS/-PKS							
<i>Single Channel w/ Jet Ventilation Channel</i>		CLM-110-HU0080 (x1), CLM-110-HU0136 (x2), and KS Connector Set 13991CK		DSD-110-HU0163, DSD-110-HU0113, and KS Connector Set 13991CK		2-8-640HAN	
Gastrosopes	13807 -NKS/-PKS						
Colonoscopes	13925 -NKS/-PKS						
<i>Dual Channels w/ Jet Ventilation Channel</i>		CLM-110-HU0080 (x1), CLM-110-HU0136 (x2), CLM-110-HU0132 (x1), and KS Connector Set 13991CK	CLM-110-HU0080 (x3), CLM-110-HU0136 (x1), CLM-110-HU0132 (x1), and KS Connector Set 13991CK	DSD-110-HU0166 and KS Connector Set 13991CK		2-8-740HAN	
Gastrosopes	13806 -NKS/-PKS						
Colonoscopes	13910 -NKS/-PKS						
	13911 -NKS/-PKS						
KARL STORZ Recommended HLD Agents for Medivators AERs:		CIDEX, CIDEX OPA, or Rapicide		Rapicide PA	Rapicide	Rapicide PA	

NOTE: The Countertop Endoscope Reprocessor (CER), Dual Sided Disinfector (DSD), Single Sided Disinfector (SSD), and DSD Edge Hook-up Connectors require the KARL STORZ Washer/Disinfector Connector Set 13991CK.