

# SCOPE BUDDY™ PLUS

## Endoscope Flushing Aid | Product Re-Order Guide

### Tubing

**24-Hour Multi-Use Pump Tube (30/box) SB-110-HU1015**

Main tubing used to pull water from the sink or clean container, through the pump, hookup and endoscope channels. Also used to dose detergent to the sink. Discard daily, do not use for more than one (1) day.

**24-Hour Multi-Use DSD Hookup Adapter (30/box) SB-110-HU1010**

Allows for connection to DSD Hookups for easy transfer to DSD and SSD Series, except for FUJIFILM™ Endoscopes. Discard daily, do not use for more than one (1) day.

**24-Hour Multi-Use Hookup for OLYMPUS® Endoscopes (30/box) SB-110-HU1007**

Allows for connection to most OLYMPUS Gastrointestinal Endoscopes. Discard daily, do not use for more than one (1) day.

**24-Hour Multi-Use Hookup for PENTAX® Endoscopes (30/box) SB-110-HU1008**

Allows for connection to most PENTAX Gastrointestinal Endoscopes. Discard daily, do not use for more than one (1) day.

**24-Hour Multi-Use Hookup for FUJIFILM™ Endoscopes (30/box) SB-110-HU1009**

Allows for connection to most FUJIFILM Endoscopes. Discard daily, do not use for more than one (1) day.

**24-Hour Multi-Use Universal Hookup (30/box) SB-110-HU1016**

Allows for connection to endoscopes that do not connect to 24-Hour Multi-Use Hookups by using manufacturer's cleaning adapters. Also used to flush up to four (4) wire-guided dilators at once. Discard daily, do not use for more than one (1) day.

**24-Hour Multi-Use Auxiliary Detergent Out Tube (30/box) SB-110-HU1012**

Doses detergent into a separate clean container. Discard daily, do not use for more than one (1) day.

**Detergent In Tube (10/box) SB-110-HU1011**

Pulls detergent from the detergent jug into the SCOPE BUDDY™ PLUS Endoscope Flushing Aid. Replace every time the detergent jug is replaced.

### Miscellaneous

**Barcode Reader with Holder (1 each) 78400-834**

Includes Barcode Reader, Programming Instructions, USB cable and Barcode Reader Holder.

**Barcode Reader Bracket (1 each; 2 pieces) 47049-880 (Bracket attaches to back of unit) & 47049-878 (Bracket end piece)**

Supplied with unit to hold barcode reader. Can be mounted to either the left or right side of the unit as desired.

**Temperature Probe (1 each) 37590-801**

The system includes a reusable temperature probe with two sensors to monitor the temperature of the detergent water. Wipe down daily with disinfect wipe.

**Wall Mounting Bracket (1 each) 72233-214**

Included with the system and is used to mount SCOPE BUDDY PLUS Endoscope Flushing Aid on the wall.

**Graduated Cylinder (1 each) 75227-679**

Used in measuring detergent for daily Dosing Verification and Calibration.

**OLYMPUS Air/Water Blanking Cap (1 each) SB-110-HU1017**

Used to block off OLYMPUS Endoscope Air/Water Ports.

**Barcode Tags (Numbered 100-199) 67199-157; (Numbered 200-299) 67199-158; (Numbered 300-399) 67199-541**

Adhesive barcode tags: includes duplicates of each number in a resealable vial.

### Detergent and Wipes

**INTERCEPT™ Detergent (4 - one gallon bottles/case) ML02-0106**

Endoscope detergent for pre-cleaning and provides superior removal of biological and organic soils. Effectively penetrates bacterial biofilms found on endoscope surfaces.

**INTERCEPT™ Wipes (4 containers/case, each container has 100 wipes) ML02-0107**

INTERCEPT Detergent towelettes are a low cost way to wipe down contaminated endoscopes and reduce biofilm.

**MEDICLEAN EZ™ Enzymatic Detergent (4 - one-gallon bottles/case) ML02-0131**

Endoscope enzymatic detergent for pre-cleaning with advanced proteolytic cleaning action.

www.medivators.com • 1-800-444-4729

SCOPE BUDDY™, INTERCEPT™ and MEDICLEAN EZ™ are trademarks of Medivators Inc.

OLYMPUS® is a registered trademark of Olympus Corporation.

PENTAX® is a registered trademark of Hoya Corporation.

FUJIFILM™ is a trademark of FUJIFILM Corporation.

© 2019 Medivators Inc. All rights reserved. 50098-1493-EN Rev C

